




**INVENT BIO-MED Pvt. Ltd.**  
 Invention Enhances Life...  
 An ISO 9001-2008, ISO 13485-2003 Certified Company




# Introducer

## ADVANCE GUIDEWIRE INTRODUCER



Bevel Tip  
(BV)



Back cut Tip  
(BC)

- Transparent hub for visualizing the guidewire.
- Available with echogenic tip.
- Smooth hub to cannula transition for easy introduction.
- Sharp cannula provides quick access.
- Bevel /Back cut tip guidewire introducer protected by plastic PP Tube.
- SS material AISI304.



### ordering information

Cat Code	Size	Length
110501 (BV/BC)	18 G	7 cm
110501 (WBC)	18 G	7 cm
110502	20 G	4 cm
110503	21 G	4 cm

**Invent Bio-med Pvt. Ltd. (Head Office)**  
 517, Laxmi Plaza, New Link Road, Andheri (W)  
 Mumbai - 400 053 INDIA.  
**Tel.:** +91-22- 2632- 4234 **Fax:** +91-22- 2632- 4250  
**Email:** sales@inventbio-med.com