



**INVENT BIO-MED Pvt. Ltd.**  
Invention Enhances Life...  
An ISO 9001-2008, ISO 13485-2003 Certified Company



# Galaxy

## Inflation Device



- Quick Latching mechanism
- 30 ATM manometer
- 20cc syringe with double lip seal
- Nylon reinforced hose
- Stable ergonomic design
- Clear body with draft - free bore

### Ordering Information

Cat Code	Description
111205	30 ATM Galaxy Inflation Device
111206	30 ATM Galaxy Inflation Device Kit
<b>Kit*</b> : Hemostasis Y-Connector + Torquer (P/W) + Blunt Introducer	

**note: \*customized kit available on request**

**Invent Bio-med Pvt. Ltd. (Head Office)**  
517, Laxmi Plaza, New Link Road, Andheri (W)  
Mumbai - 400 053 INDIA.  
**Tel.:** +91-22- 2632- 4234 **Fax:** +91-22- 2632- 4250  
**Email:** sales@inventbio-med.com