




INVENT BIO-MED Pvt. Ltd.
Invention Enhances Life...
An ISO 9001-2008, ISO 13485-2003 Certified Company



Introducer

SHEATH



- Radio opaque sheath with proprietary homeostasis valve and side apart with attached 3-way stopcock side arm.
- Radio opaque atraumatic vessel dilator.
- Smooth atraumatic arterial sheath & dilator transitions locking hub design.

Ordering Information

Cat Code	French Size	Length*	Guidewire	GW Introducer
111101	4F	10cm	0.038"	18G
111102	5F	11cm	0.038"	18G
111103	6F	11cm	0.038"	18G
111104	7F	11cm	0.038"	18G
111105	8F	11cm	0.038"	18G
111106	9F	11cm	0.038"	18G
111107	5F	23cm	0.038"	18G
111108	6F	23cm	0.038"	18G
111109	7F	23cm	0.038"	18G
111110	8F	23cm	0.038"	18G
111111	9F	23cm	0.038"	18G

note: *can be tailored as per customer requirement

Invent Bio-med Pvt. Ltd. (Head Office)
517, Laxmi Plaza, New Link Road, Andheri (W)
Mumbai - 400 053 INDIA.
Tel.: +91-22- 2632- 4234 **Fax:** +91-22- 2632- 4250
Email: sales@inventbio-med.com