




INVENT BIO-MED Pvt. Ltd.
Invention Enhances Life...
An ISO 9001-2008, ISO 13485-2003 Certified Company



NOVA

Inflation Device

- 30 atm gauge reliable in any procedure.
- 90 degree rotating manometer allows convenience in handling.
- Coded figures from nominal to burst pressure.
- Unique mechanism for rapid pressure increase helps to increase efficiency for fast procedures.
- 20cc graduated compact syringe.
- Grip design for better handling.
- Clear material for easy air bubble detection and removal.
- Ergonomic design for optimal handling.
- Progressive resistance while inflation gives pressure feedback.



Ordering Information

Cat Code	Description
111207	30 ATM NOVA Inflation Device
111208	30 ATM NOVA Inflation Device Kit
Kit* : Hemostasis Y-Connector + Torquer (P/W) + Blunt Introducer	

note: *customized kit available on request

Invent Bio-med Pvt. Ltd. (Head Office)
517, Laxmi Plaza, New Link Road, Andheri (W)
Mumbai - 400 053 INDIA.
Tel.: +91-22- 2632- 4234 **Fax:** +91-22- 2632- 4250
Email: sales@inventbio-med.com