




**INVENT BIO-MED** Pvt. Ltd.  
Invention Enhances Life...  
An ISO 9001-2008, ISO 13485-2003 Certified Company



# Custom Kit

## PTCA & ANGIOGRAPHY



- Large Bore 9F 'Y' Connector, Insertion Tool and Torquer to facilitate the PTCA / STENT procedure.
- High Pressure 1200 PSI Angiodye Connector for easy flow of contrast.
- Control Device 12cc, Transfer Set.
- PM. Line Connector M/F, M/M.

### Ordering Information

Cat Code	Description
110902	PTCA KIT 5* items
2 port- 3 port <b>Manifolds</b> , Control Device, <b>Y-Connector</b> , Pressure Monitoring Line & Insertion tool	

**Invent Bio-med Pvt. Ltd. (Head Office)**  
517, Laxmi Plaza, New Link Road, Andheri (W)  
Mumbai - 400 053 INDIA.  
**Tel.:** +91-22- 2632- 4234 **Fax:** +91-22- 2632- 4250  
**Email:** sales@inventbio-med.com