



INVENT BIO-MED Pvt. Ltd.
Invention Enhances Life...
An ISO 9001-2008, ISO 13485-2003 Certified Company



Push Click



HEMOSTASIS Y-CONNECTOR

- For easy, fast and accurate single handed intra vascular device insertion.
- Perfect insertion device respecting the integrity of the coated surface.
- Combination of a hemostatic valve and an integrated push-pull valve opener.
- Clean and safe procedure with a minimized backflow thus limiting blood loss.



Ordering Information

Cat Code	Description
110306	Hemostasis Push Click Y- Connector
110307	Hemostasis Push Click Y- Connector Kit with side arm stopcock
Kit* = Hemostasis Push Click Y- Connector + Torquer (P/W) + Blunt Introducer	
110904	Hemostasis Push Click Y- Connector Kit without side arm stopcock
Kit* = Hemostasis Y- Connector + Torquer (P/W)+ Blunt Introducer	

note: *customized kit available on request

Invent Bio-med Pvt. Ltd. (Head Office)
517, Laxmi Plaza, New Link Road, Andheri (W)
Mumbai - 400 053 INDIA.
Tel.: +91-22- 2632- 4234 **Fax:** +91-22- 2632- 4250
Email: sales@inventbio-med.com